

**Tabouret de douche Oasis**



The Oasis shower stool allows you to sit comfortably in a shower tray, even from reduced size. The stool has armrests, very convenient for helping the elderly to stand up or

reduced mobility.

The seat is removable for easy cleaning. The large width between the armrests (52 cm) allows the seat

strong people. Ultra-stable and resistant stool thanks to

to its reinforced structure.

Floor dimensions: width. 48 x prof. 47 cm.

Dimensions of

the seat: width. 42 x prof. 28 cm. Height

adjustable seat

from 48 to 63 cm.

Maximum authorized weight 154 kg.