Days Chrome Plated Wheelchairs

Self-Propelled & Attendant Wheelchair



41cm (16")

21823FBN

23823FBN

Folding Back

Designed for strength and durability, this quality range of wheelchairs are manufactured from heavy gauge steel and have an attractive scratch-resistant chrome finish. The embossed vinyl upholstery is hard-wearing and easy to wipe clean. Puncture-proof tyres ensure low maintenance, yet deliver a smooth ride over the majority of surfaces.

Supplied as standard with padded full length arms for comfort and support, desk length arms can be added to enable easier access to tables and worktops. The height adjustable, swing-away and removable footrests come complete with a calf strap and can be easily replaced with multi-adjustable elevating leg rests.


- · Available in a choice of two seat widths
- · Detachable swing-away armrests and footrests
- · Folds easily for storage and transport
- · Moulded armrests for extra comfort
- · Chrome plated steel frame for durability
- Stainless steel skirt guards
- Wipe clean padded seat fabric
- · Available with fixed or folding backrest options
- · Optional elevating leg rests, desk arms and adjustable I.V. pole

Seat Width	Narrow: 41cm (16"), Standard: 46cm (18")
Seat Depth	41cm (16")
Seat Height (Front)	51cm (20")
Seat Height (Rear)	47cm (181/2")
Backrest Height Above Seat	41cm (16")
Footrest Length	38 - 53cm (15 - 21")
Armrest Height Above Seat	25cm (9¾")

	Self-Propelled
External Width (Open)	62cm (211/2"), 67cm (261/4")
External Width (Closed)	30cm (11¾")
External Depth (w/o Footrests)	80cm (311/2")
External Depth (w Footrests)	105cm (41¼")
Chair Weight	20kg (44lb)
Chair Weight (w/o Footrests)	18kg (39.6lb)

Elevating Leg Rests

RP36 RP36N





Cast Support RP70L Left RP70R Right Desk Arms RP54

AYear Main Frame