

0

HK 54 Cosy Shoulder & neck heat pad


engineering



Machine washable up to 30°C

Cosy and soft heat pad

Formfitted form for shoulder and neck	
Breathable, cosy and skin-friendly	
Cosy wearing comfort - with magnet closure	
ndividually adjustable by press studs	
Electronic temperature control	
3 illuminated temperature settings	
Automatic switch off after approx. 90 minutes	
Rapid heat up	
3SS (Beurer Safety System)	
Detachable control unit	
Deko-Tex	
he textiles used for this device meet the stringent hu	uman
ecological requirements of Oeko-Tex Standard 100,	
as verified by Hohenstein Research Institute.	
1ains operation:	
220-240V; 50-60 Hz; 100 W	
Product measurements (LxB): 56 x 52 cm	
Product weight: approx. 770 g	
3 year guarantee	
Sales unit: 4 / Shipping carton: –	
AN no.: 4211125 21616 0	
tem no.: 216.16	



216.16_HK54-cosy_2019-06-14_01_DS_EN_BEU Subject to errors and changes