

Product Summary

Product Name: G4 Integrated Bariatric Hybrid Mattress system

Product Code: 8248

<u>Design Platform:</u> Integrated Hybrid Mattress system (foam-in-cell)

Specification summary:

Level of risk - Very High

Mattress Dimensions – 2,000 x 1200 x 150mm

Mattress net weight – 22kg

Pump unit dimensions – 265 x 175 x 90 mm

Pump net weight – 1.7kg

User weight threshold - Up to 60 stone / 380 kg (no lower specification limit)

Functionality: Cycle completion – 12 minutes

Electronic pressure sensing control

Voltage – 220-240V, 50 Hz

Power Consumption - 10W

2 in 1 cell inflation

Technical support:

Warranty - 2 years

Warranty extension – yes, when manufacturer-serviced annually

Annual service – yes, carried out by UK-based manufacturer

Spare parts availability – readily available, UK-based manufacturer.

Replacement covers – Choice of HF-welded or stitched, full or section zips

Product care/maintenance:

Infection control – Integrated cable management

- Integrated air hose cover which is wipeable and machine-washable
- 1,000 ppm diluted or 0.1% Hypochlorite solution (or to be used in conjunction with local infection control guidelines)

Replacements - Individually replaceable top and bottom covers

- ALL air-pump key assemblies available as spare parts to order

Features & Benefits:

- Combination of foam and air used by the integrated hybrid mattress design allows for broader clinical application up to very high risk level
- Design for easy profile when used with profiling bed frames
- Integrated carry handles for easy transportation and storage
- Integrated visual and audible alarms when faults detected
- Rapid deflation via integrated CPR
- Low and high pressure settings
- Unique, patented air-cell design which prevents air cell twisting