

The Alliana Seat Cushion Affordable pressure care



- ✓ Soft PVC Air Cells
- ✓ Adjustable comfort
- All seams fully welded
- Belgian two way stretch covers
- Zip closure for easy access & cleaning

Specifications

Alliana Seat Cushion

Dimensions Standa	ard L 48 x W 44 x H 6cm
Cell Type	56 x PVC Soft Single Mould Cells
Top Cover	Dahlia Bi-Elastic 2 way stretch - welded
Base Standard	Dahlia Bi-Elastic 2 way stretch - welded

Order Code: ALLI-STC-44x48x6

Inflation	ŀ	Hand pump and valve
Infection Control	All seams fully welde	d - waterfall over zip
Base Attachment Standard		Rear Velcro Straps
Warranty	2 years Soft Goods	ARTG: 175810

Unit C3, 15 Narabang Way Belrose, Sydney NSW 2085 t: 02 9844 5456 f: 02 9844 5445 e: info@patientsupportsystems.com

www.patientsupportsystems.com

The Alliana Air Alternating Cushion BlueTooth 3 in 1 Therapy Control





Unit C3, 15 Narabang Way Belrose, Sydney NSW 2085 **1800 432 247** e: info@patientsupportsystems.com www.patientsupportsystems.com

