

TL 30

beurer
wellbeing

TL 30 Daylight therapy lamp

 **german
engineering**



Simulation of daylight:
Light intensity of approx. 10,000 lux (distance 10 cm*)



Energy-saving LED technology



Medical device

For use in case of light deprivation symptoms:
e.g. low mood, lack of energy and drive

Compact size

Ideal for your desk or on the way - Practical travel companion

Exceptionally even illumination

Light intensity: approx. 10,000 lux (distance 10 cm*)

Flicker-free & UV-free

Convenient one-button operation

Versatile stand

- horizontal or vertical applicable
- continuously adjustable

Incl. mains adapter, cable length approx. 148 cm

Incl. practical storage bag

Medical device

Colour rendering index CRI: > 80

Colour rendering index (CRI) describes the realistic reproduction and perception of colours under artificial light compared to natural light. CRI is a quality feature for artificial light sources (sunlight CRI = 100).

Color temperature: 6,000 - 6,500 K

Illumination surface: approx. 12 x 20 cm

Product measurements: approx. 15.6 x 2.6 x 23.6 cm


Product weight: approx. 320 g

3 year guarantee

Sales unit: 4 / Shipping carton: -

EAN-Nr.: 4211125608057

Art.-Nr.: 608.05



| * Lux | Abstand / Distance / Distance Distanza / Distanza |
|--------|--|
| 10.000 | ca./approx./env./aprox./ ca. 10 cm |
| 5.000 | ca./approx./env./aprox./ ca. 20 cm |
| 2.500 | ca./approx./env./aprox./ ca. 30 cm |

** without harmful UV radiation

608.05_TL30_2020-09-22_02_DS_EN_BEU Subject to errors and changes