

WOBURN ULTRA-LOW BED





Product Code: HLB799

Description: • A fully functioning four- section Ultra-Low profiling bed.

• Platform can be lowered to 67mm off the floor for patient safety.

• Electronically setting of height, back rest or knee brake is simple, using the lockable handset.

• The Woburn Ultra-Low has a Trendelenburg function as standard.

• Optional bed levels, side rails, lighting and accessories available

Safe Working Load: 185 kg (29 Stone)

Overall Dimension: 2130 x 900mm

Adjustable Platform Height: 67mm -640mm

Mattress Dimensions: 2000 x 900mm

Angle of Back Rest: Up to 70 degrees

Angle of Knee Brake Up to 30 degrees

Bed Weight: 81 kg

Key Information:Manual Leg Rest Adjustment

Wooden Head and Foot Boards- Wrapped for decontamination

Low Voltage- 29V; 9V Battery Backup

Easy Cleaning- Complies to MHRA Guidelines

Comprehensive User Manual

CE Marked

Subject to the terms of Harvest Healthcare Guarantee and Warranties the Ultra Low Bed is guaranteed for five years from date of purchase.

Supplier Contact Details:

t: +44(0)1709 377172 e: enquiries@harvesthealthcare.co.uk w: www.harvesthealthcare.co.uk Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60

