

Product Data Sheet

27/3/2018

Mobile Infusion Stands - Stainless Steel

FEATURES

- Meets appendix A (Stability Requirements) of BS 3619:1976 Mobile Infusion Stands
- Seamless welded 304 grade stainless steel construction
- Height adjustable upright using a handwheel (1270 - 2000mm)
- Available with two or four pigtail hooks
- 9kg plastic coated weighted base available for additional stability
- Maximum recommended load:
- Hooks - 4kg per stand plus
- Fixed Upright
- Non-Weighted Base - 2kg per stand
- Weighted Base - 4kg per stand
- 50mm polypropylene/thermoplastic rubber twin wheel castors (1 x anti-static)

STANDARDS

- Meets appendix A (Stability Requirements) of BS 3619:1976 Mobile Infusion Stands

DIMENSIONS (w x d x h)

- Overall Dimensions - 650 x 650 x 1270-2000mm

STANDARD COLOURS

- Base/Upright - Stainless Steel
- Weighted Base - Grey

PRODUCT LINK

http://www.bristolmaid.com/product.asp?CAT_ID=178&P_ID=386



DPWS/1



DPWS/1/W9/4H

Item	Description	Weight
DPWS/1	Two Hook	3.8kg
DPWS/1/4H	Four Hook	3.9kg
DPWS/1/W9	Two Hook, Weighted Base	12.8kg
DPWS/1/W9/4H	Four Hook, Weighted Base	12.9kg

For information on recycling & cleaning instructions please refer to www.bristolmaid.com or contact our technical services department on 01258 451338

Blandford Heights, Blandford Forum, Dorset DT11 7TG, UK

Tel: 01258 451338

Fax: 01258 455056

Email: sales@bristolmaid.com

www.bristolmaid.com



Bristol MaidTM
Hospital Metalcraft Limited 