



HAND CONTROL T4



FOOT CONTROL





EXPLANATION OF OTHER PICTOGRAMS

Pictogram	Denomination	Description
\triangle	Warning	This sign calls for your special attention.
C€	CE-mark	The mark states that CE-marking for the system has been elaborated.
F	IMPORTANT	Please note special information.
	Electronic waste	The symbol indicates that the system upon terminated use should be delivered to the nearest recycling and disposal according to ruling laws and regulations.
	Level of sealing	Level of sealing IPX0 => No protection against water destroying the device. Not suitable for outdoor use.

THE THREE POSITIONS

For step-by-step instructions – please see chapter 5 in the manual

1: Working position (during turning)





2: Folded position



(Out of bed position)

3: Bed rail position



(during turning)