



## IH 60 Nebuliser



Nebulisation of liquid medicine using compressed air technology



Rechargeable, with lithium-ion battery



Very quiet, small and light

For the treatment of the upper and lower respiratory tract For use in colds, asthma, respiratory disease, etc.

- for fast and effective inhalation

High proportion of respirable particles

High nebulisation capacity (approx. 0.25 ml/min) - short inhalation time

Accessories:

Nebuliser bottle, mouthpiece, nosepiece, adult mask, child mask, compressed air tube, mains adapter with micro-USB cable, filter, storage bag and carrying bag

Can be desinfected

Operating pressure: approx. 0.25-0.5 bar

Medical device

Li-ion battery: 2570 mAh; 3.7 V

Mains operation:

Input: 100-240 V~; 50-60 Hz; 0.5 A

Output: 5 V === 2 A

Product measurements (LxWxH): 110 x 62 x 47 mm

Product weight: approx. 222 g

5 year guarantee

Sales unit: 3 / Shipping carton: 4 x 3

EAN no.: 4211125 60206 2

Item no.: 602.06



german